

# CENTURY

SERIES Preferred by medical facilities worldwide.



Optional adjustable contoured backrest with black plastic shroud

16" wide, 4" thick poly-foam seat cushion


Pneumatic or precision spinlift height adjustment

Optional adjustable chrome foot ring

Chrome hooded dual wheel casters

Five-leg, 23" brushed aluminum accent base

height range



B - backrest  
 V - seamless/vacuum formed  
 LCD - locking casters  
 G - glides  
 FR - foot ring

\*Height range 16.75" – 23.5"