

VALUE SERIES MAT PLATFORMS

	Width	Length	Height
269-47	4'	7'	18"
269-57	5'	7'	18"
269-68	6'	8'	18"

All Laminate Upholstered Top Mat Platform

- 2" foam padding
- All Gray laminate construction
- 3 full panel legs
- All hardwood Plywood frame
- Easy-clean surfaces
- Solid, bolt-through construction
- Frame attached with steel angles
- 500 lbs. load capacity



	Width	Length	Height
260-47	4'	7'	18"
260-57	5'	7'	18"
260-68	6'	8'	18"

Upholstered Top Mat Platform

- 2" foam padding
- 5 triple-bolted, solid hardwood legs
- All hardwood Plywood frame
- Solid, bolt-through construction
- Frame attached with steel angles
- 500 lbs. load capacity

	Width	Length	Height
262-47	4'	7'	18"
262-57	5'	7'	18"
262-68	6'	8'	18"

Mat Platform with removable Mat

- Mat attached with Velcro® strips to the rimless subtop
- 2" mat, in any of our vinyl colors, included
- 5 triple-bolted, solid hardwood legs
- All hardwood Plywood frame
- Solid, bolt-through construction
- Frame attached with steel angles
- 500 lbs. load capacity



Mat attached with Velcro® strips

CLASSIC WOOD MAT PLATFORMS

FEATURES

- 18" standard height with 21" height option
- 16 standard colors of upholstery Perma-Block-3 Protection
- Available in Red Oak as an option (see options below)
- Bariatric load capacities (see specific models for individual load capacities)
- 4 models to choose from


	Width	Length	Height
240-47	4'	7'	18"
240-57	5'	7'	18"
240-68	6'	8'	18"

Upholstered Mat Platform

- Features all hardwood construction
- 700 lbs. load capacity

Options (Classic Series Only)

-042	Red Oak
-041	21" Height



CLASSIC WOOD MAT PLATFORMS



	Width	Length	Height	Backrest Dimensions
241-47	4'	7'	18"	4' x 30"
241-57	5'	7'	18"	5' x 30"
241-68	6'	8'	18"	3' x 30" each side

Upholstered Mat Platform with adjustable Backrest

- Features 6-position, all-steel, adjustable backrest
- 6' x 8' table has a split backrest
- 500 lbs. load capacity



Mat attached with Velcro® strips



	Width	Length	Height
242-47	4'	7'	18"
242-57	5'	7'	18"
242-68	6'	8'	18"

Mat Platform with Removable Mat

- Mat attaches with Velcro® strips to the rimless subtop
- 2" mat, in any of our vinyl colors, included
- 700 lbs. load capacity
- Mat included



	Width	Length	Height
243-47	4'	7'	18"
243-57	5'	7'	18"

Space Saving, Folding Mat Platform

- Features full back wall that is floor supported
- Sturdy 3 front-leg design
- Projects only 8" from wall when folded
- 500 lbs. load capacity



ALPHA SERIES MAT PLATFORMS

FEATURES

- All-welded, heavy duty steel frame finished in Gray
- Heavy-duty 1 1/2" x 3" steel legs
- Plywood top
- Solid six-leg support
- 18" fixed height—standard
- 1,000 lbs. weight capacity



	Width	Length	Height
230-47	4'	7'	18"
230-57	5'	7'	18"
230-68	6'	8'	18"

Upholstered Top Mat Platform

- 2" of firm foam padding covered with heavy weight vinyl with Perma-Block top coating
- Custom edge bumper guard helps to prevent top from tearing



	Width	Length	Height
231-47	4'	7'	18"
231-57	5'	7'	18"
231-68	6'	8'	18"

Upholstered Top Mat Platform with Adjustable Backrest

- 2" of firm foam padding covered with heavy weight vinyl with Perma-Block top coating
- Custom edge bumper guard helps to prevent top from tearing
- 6-position, all-steel, adjustable backrest
- 6' x 8' model has split 30" high backrest



	Width	Length	Height
232-47	4'	7'	18"
232-57	5'	7'	18"
232-68	6'	8'	18"

Mat Platform with Removable Top

- 2" mat attaches with Velcro® strips to the rimless subtop
- Mat is included



Mat attached with Velcro® strips

HI-LO MAT PLATFORMS



	Width	Length	Height Range
251-47	4'	7'	23"-35"
251-57	5'	7'	23"-35"
251-68	6'	8'	23"-35"

Hand-Crank Hydraulic Hi-Lo Mat Platform

- All-welded steel frame with laminated finish sides
- Stable 4-leg lift system with levelers
- 2" thick firm padding
- Custom bumper strip protects top's upholstery
- Easy to operate hand-crank height adjustment
- Hand crank easily folds for storage
- 750 lbs. lift capacity

	Width	Length	Height Range
253-47	4'	7'	22"-34"
253-57	5'	7'	22"-34"
253-68	6'	8'	22"-34"

Hi-Lo Mat Platform

- All-welded steel frame
- Hands-free, foot peddle control
- Dependable single motor drive
- 110 volts, 60 hz., 2.7 amps. with overload protection
- Fully enclosed drive unit for safety
- Custom bumper strip protects top's upholstery
- 600 lbs. lift capacity



Mat attached with Velcro® strips

	Width	Length	Height Range
252-47	4'	7'	22"-34"
252-57	5'	7'	22"-34"
252-68	6'	8'	22"-34"

Hi-Lo Mat Platform with Removable Mat

- All-welded steel frame
- Hands-free, foot peddle control
- Dependable single motor drive
- 110 volts, 60 hz., 2.7 amps. with overload protection
- Fully enclosed drive unit for safety
- 600 lbs. lift capacity
- Removable 2" mat held in place with Velcro® strips
- Mat included