



Aneroid Blood Pressure Kits



Technologically advanced yet easy to use blood pressure monitors for consumers (UA-100, left) and professionals (UA-200, right)

The UA-100 and UA-200 Aneroid Blood Pressure Kits provide consumers and health care professionals with reliable monitoring products at home or in a clinical setting. These aneroid blood pressure kits can assist in the management of hypertension, improve patient compliance to treatment of high blood pressure, and be used as tools in a preventative health management program.

- Nylon cuff for superior durability and appearance
- Latex bulb with standard air release valve
- Includes stethoscope (Model UA-100 only)
- Two year limited warranty for manometer accuracy

- Professional accuracy
- Includes zippered nylon storage case
- Calibration screwdriver included



UA-100



UA-200

For questions about **LifeSource** monitors, call our toll free **LifeSource** HealthLine at 1-888-726-9966 and a specially trained representative will assist you, or visit us at www.LifeSourceOnline.com.



Aneroid Blood Pressure Kits

A&D Medical offers two aneroid blood pressure kits. The Home Blood Pressure Kit with Attached Stethoscope (Model UA-100) is designed to give consumers a way to track their readings at home. The Professional Sphygmomanometer (Model UA-200) is designed for health professionals to measure a patient's blood pressure.

	PART NUMBER	UPC NUMBER	CUFF SIZE
UNITED STATES AND CANADA	UA-100	0-93764-00100-0	Adult 11.5" - 16" (29-41cm)
	UA-200	0-93764-00200-7	Adult 11.5" - 16" (29-41cm)

Home Blood Pressure Kit with Attached Stethoscope Model UA-100



- No-pin-stop manometer with easy to calibrate feature
- Comes with attached stethoscope and D-ring cuff
- Latex 2-tube bladder
- Includes trilingual instruction manual
- Zippered nylon storage case and calibration screwdriver included

Professional Sphygmomanometer Model UA-200



- Arrow indicator on cuff ensures proper cuff placement
- Comes with a standard adult cuff size. Four additional cuff sizes are available (see below)
- Zippered nylon storage case and calibration screwdriver included

Accessories

PART #	DESCRIPTION
11-001	Inflator Bulb: UA-100, UA-200
11-002	Exhaust Valve: UA-100, UA-200
11-003	End Valve: UA-100, UA-200
11-004	Gauge: UA-100, UA-200
11-005	Zippered Case: UA-100, UA-200
UA-700CH	Cuff Hose: UA-100, UA-200

CUFFS AND BLADDERS FOR UA-200 PROFESSIONAL SPHYGMOMANOMETERS

PART #	DESCRIPTION
11-008	Infant Cuff w/ bladder: 5.5" - 7.5" (14-19cm)
11-009	Infant bladder
11-010	Child Cuff w/ bladder: 7.5" - 10.5" (19-27cm)
11-011	Child bladder
11-006	Adult Cuff w/ bladder: 11.5" - 16" (29-41cm)
11-007	Adult bladder
11-012	Large Cuff w/ bladder: 13.5" - 20" (35-51cm)
11-013	Large bladder
11-014	Thigh Cuff w/ bladder: 16" - 23" (41-58cm)
11-015	Thigh bladder

Packaging and Shipping Specifications

Selling Unit (Box)	L x W x H (inches)	WT (lbs.)	Individual Item Cube (cubic feet)
UA-100	4 3/4 x 2 3/16 x 7 1/2	1.30	0.06
UA-200	4 3/4 x 2 3/16 x 7 1/2	1.00	0.06
Unit Case (20 units)	L x W x H (inches)	WT (lbs.)	Case Cube (cubic feet)
UA-100	17 x 11 x 16	30.10	1.73
UA-200	17 x 11 x 16	22.90	1.73



1555 McCandless Drive
Milpitas, CA 95035
www.LifeSourceOnline.com

LifeSource
Your source for a lifetime of health



LifeSource is a trademark of A&D Medical. ©2002 A&D Medical. All rights reserved.



Ask your A&D Sales Representative about other **LifeSource** products.

- Blood Pressure Monitors
- Heart Rate Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes

For questions about **LifeSource** monitors, call our toll-free **LifeSource** HealthLine at 1-888-726-9966 or visit our web site at www.LifeSourceOnline.com

LIT:UA-100 06/02