



UA-774 Dual Memory Auto-Inflate



- Ideal for two or more people
- 60 memory capacity
- Provides average blood pressure readings
- AC Adapter included
- Two cuffs included (Medium and Large)
- Latex free



A complete blood pressure kit for the entire family – a great value!

Keeping track of blood pressure readings for the entire family can now be done with just one monitor. The UA-774 Dual Memory Auto-Inflate blood pressure monitor provides individual and average blood pressure readings for two. A guest mode is also included for others to obtain their readings. Included with this monitor are two cuffs (medium and large) and an AC Adapter.

This unique monitor is equipped with two start buttons - one for each user - and each button can store up to 30 readings for a total of 60 measurements. The monitor also features A&D's unique Irregular Heartbeat Feature which provides a blood pressure and pulse rate readings even when an irregular heartbeat occurs, making this an ideal kit for the entire family.

- Modern, illuminated display
- Offers great value to customer
- Irregular Heartbeat feature
- Large-print trilingual instruction guide with illustrations and quick start card in English, French and Spanish
- Four week blood pressure logbook



For questions about *LifeSource* monitors, call our toll free *LifeSource* Health Line at 1-800-726-3364 and a specially trained representative will assist you, or visit us at www.LifeSourceOnline.com.

UA-774 Dual Memory Auto-Inflate

To obtain an accurate reading, it is important to use the correct cuff size. A cuff that is too small will result in a higher than expected reading and a cuff that is too large will result in a lower than expected reading. The UA-774 Dual Memory Auto-Inflate unit includes two different cuff sizes to accommodate your customer's needs.

	PART NUMBER	UPC NUMBER	CUFF SIZES
UNITED STATES AND CANADA	UA-774AC	0-93764-60108-8	Medium 7.5 - 12.2" (19 - 31 cm)
			Large 11.8 - 17.7" (30 - 45 cm)

Each unit comes with a trilingual instruction manual, quick start card, four week blood pressure logbook and warranty card.

Accessories

PART #	DESCRIPTION
UA-279	Small Blood Pressure Monitor Cuff: 5.1 - 7.9" (13 - 20 cm)
UA-280	Medium Blood Pressure Monitor Cuff: 7.5 - 12.2" (19 - 31 cm)
UA-281	Large Blood Pressure Monitor Cuff: 11.8 - 17.7" (30 - 45 cm)
UA-380	Medium <i>EasyCuff</i> ™ Feature: 7.5 - 14.2" (19 - 36 cm)
UA-700CH	Cuff Hose
TB:181	AC Power Adapter

Technical Specifications

Dimensions	4.3" (l) x 5.8" (w) x 2.5" (h) (110 mm x 147 mm x 64 mm)
Weight	10.6 oz. (300 g) without batteries
Display	Digital character height: 16 mm (blood pressure), 10 mm (pulse rate) Pressure and pulse displayed simultaneously
Memory	60 readings (30 for each start button)
Pressurization	Automatic, using micropump
Deflation	Automatic exhaust
Type	Oscillometric
Accuracy	Pressure: ±3 mmHg or ±2%, whichever is greater; Pulse: ±5%
Pressure Range	20 mmHg to 280 mmHg; Pulse: 40 to 200 pulses/minute
Power Source	Four type AA alkaline batteries (1.5 volt) or AC Power Adapter (TB:181). Batteries not included.
Battery Life	Approx. six months with one daily measurement
Operating Environ.	50° to 104° F (10° to 40° C)
Storage Environ.	14° to 140° F (-10° to 60° C)

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method within the limits prescribed by the American National Standards Institute (ANSI/AAMI SP10) for electronic sphygmomanometers. Manufactured by an ISO 9000 and ISO 13485 registered company. Engineered in Japan. Made in China.

Packaging and Shipping Specifications

Selling Unit (Box)	L x W x H (inches)	WT (lb.)	Individual Item Cube (cubic feet)
UA-774	5 1/8 x 6 3/4 x 8 1/8	2.76	0.16
Unit Case (10 units)	L x W x H (inches)	WT (lb.)	Case Cube (cubic feet)
UA-774	26 x 13 3/4 x 8 1/2	29.78	1.71

Ask your A&D Sales Representative about point-of-purchase floor and counter displays available for this blood pressure monitor unit.



1555 McCandless Drive
Milpitas, CA 95035
www.LifeSourceOnline.com

LifeSource™
Your source for a lifetime of health



LifeSource and EasyCuff are trademarks of A&D Medical. ©2003 A&D Medical. All rights reserved.



Ask your A&D Sales Representative about other *LifeSource* products.

- Blood Pressure Monitors
- Heart Rate Monitors
- Nebulizers
- Digital Thermometers
- Personal Scales
- Stethoscopes

For questions about *LifeSource* monitors, call our toll-free *LifeSource* Health Line at 1-800-726-3364 or visit our web site at www.LifeSourceOnline.com

LIFLS774 11/03