

DUAL MEMORY WRIST UB-512 Blood Pressure Monitor

Fast. Easy. Accurate.



An advanced wrist monitor for two or more people.

A&D Medical has developed state-of-the-art measurement technology to monitor blood pressure at the wrist for two or more people. This wrist monitor is smaller, lighter and more compact than other wrist monitors and is ideal for travel. Unlike other monitors that measure blood pressure during deflation of the cuff, the Dual Memory Wrist monitor measures blood pressure during cuff *inflation* for a faster, more comfortable reading. Each START button stores 30 blood pressure measurements for a total of 60 readings.

The average reading feature automatically calculates the “average” of the total readings stored in memory. The Dual Memory Wrist monitor also offers two innovative features, the Pressure Rating Indicator™ and Irregular Heartbeat feature. The Pressure Rating Indicator provides a snapshot of your blood pressure classification based on your measurements. In addition, this monitor alerts the user of the presence of an Irregular Heartbeat and provides an accurate blood pressure and pulse rate measurements even when an Irregular Heartbeat occurs.

- Includes two “AAA” batteries
- Large-print instruction guide with illustrations and quick start card in English, French and Spanish
- Four week blood pressure logbook included
- Carrying case included

For questions about LifeSource products, call us toll free at 1-800-726-3364 and a specially trained representative will assist you, or visit us at www.LifeSourceOnline.com.

- Professional accuracy
- For two or more people
- 60 reading memory
- Pressure Rating Indicator™
- Irregular Heartbeat feature
- Displays average readings
- Fast measurement
- Ideal for travel



Dual Memory Wrist UB-512 Blood Pressure Monitor

The Dual Memory Wrist blood pressure monitor fits easily and comfortably to one's wrist to provide an accurate reading. This compact and portable unit gives fast measurements and provides convenience to those looking for an easy way to monitor their blood pressure away from home.

	PART NUMBER	UPC NUMBER	CUFF SIZE
UNITED STATES	UB-512	0-93764-60153-8	5.3 - 8.5" (13.5 - 21.5 cm)
CANADA	UB-512CN	0-93764-60154-5	5.3 - 8.5" (13.5 - 21.5 cm)

Each unit comes with a trilingual instruction manual, quick start card, two AAA batteries, carrying case, four week blood pressure logbook and warranty card.

Technical Specifications

Dimensions	0.8 (d) x 2.5 (w) x 2.2" (h) (21 x 63 x 56 mm)
Dimensions with Case	3.5 (d) x 3.4 (w) x 2.7" (h) (89 x 86 x 69 mm)
Weight	2.9 oz. (82 g) without batteries
Display	Digital: 10/7 mm character height Pressure and pulse displayed simultaneously
Memory	60 readings
Pressurization	Automatic, using micropump
Deflation	Constant air release-valve system
Type	Oscillometric
Accuracy	Pressure: ± 3 mmHg or 2%, whichever is greater; Pulse: $\pm 5\%$
Pressure Range	20mmHg to 280 mmHg; Pulse: 40 to 200 pulses/minute
Power Source	Two type AAA alkaline batteries (1.5 volt). Batteries included.
Battery Life	Approx. 400 measurements
Operating Environ.	50° to 104°F (10° to 40°C)
Storage Environ.	14° to 140°F (-10° to 60°C)

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method within the limits prescribed by the American National Standards Institute (ANSI/AAMI SP10) for electronic sphygmomanometers. Manufactured by an ISO 9001 and ISO 13485 registered company. Engineered in Japan. Made in China.

Packaging and Shipping Specifications

Individual Unit (Box)	W x D x H (inches)	WT (lb.)	Individual Item Cube (cubic feet)
All UB-512 models	4 ^{9/16} x 3 ^{7/16} x 6 ^{1/16}	0.71	0.06
Unit Case (10 units)	W x D x H (inches)	WT (lb.)	Case Cube (cubic feet)
All UB-512 models	18 x 10 x 7	8.40	0.73

Ask your A&D Sales Representative about point-of-purchase floor and counter displays available for this blood pressure monitor unit.



1555 McCandless Drive
Milpitas, CA 95035
www.LifeSourceOnline.com

LIFE SOURCE®



Pressure Rating Indicator is a trademark of A&D Medical. Patent Pending. LifeSource is a registered trademark of A&D Medical. ©2005 A&D Medical. All rights reserved.



Ask your A&D Sales Representative about other LifeSource products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes

For questions about LifeSource products, call us toll-free at 1-800-726-3364 or visit our web site at www.LifeSourceOnline.com

LITLS512 8/05