

#1 Brand Recommended By Pharmacists and Doctors

OMRON®


HEM-712C

Automatic Blood Pressure Monitor

American Heart Association® 
Learn and LiveSM



New & Improved

-  **IntelliSense®** Technology for More Accurate, More Comfortable Readings
- **Simple, Silent, 1-Touch Automatic Operation**
- **Extra Large Digital Display**
 - Systolic and Diastolic Blood Pressure
 - Pulse
 - Date and Time
- **Advanced Diagnostics**
 - Irregular Heartbeat Detector (♡)
 - 30 Memory Storage (📄)
 - Advanced Averaging Displays Up to Last 3 Readings in Last 10 Minutes (📊)
- **Durable Contoured Arm Cuff**
 - Fits Arms 9"-13" in Circumference
(Large Cuff, 13" to 17", available separately)
 - D-ring Cuff is Easy to Apply and Fits Securely
- **Illustrated Instruction Manual Included**
- **2 Year Warranty**

A Good Sense of Health

MODEL HEM-712C

PACKAGING SPECIFICATIONS...

CARTON:	QTY:	APPROX. WEIGHT:	APPROX. DIMENSIONS: (W x D x H)
Unit: Package	1 each	1.3 lbs.	5.25" x 4.87" x 7"
Master:	12 each	18.3 lbs.	20.35" x 17.83" x 8.18"
UPC:	0 73796 07122 6		

PRODUCT SPECIFICATIONS...

Display:	Liquid Crystal Display
Measurement Range:	Pressure: 0 - 280 mmHg Pulse: 40 - 180/minute
Accuracy/Calibration:	Pressure: ± 3 mmHg or 2 % of reading Pulse: $\pm 5\%$ of reading
Inflation:	Automatic by electric pump
Deflation:	Automatic pressure release valve
Rapid Air Release:	Automatic exhaust valve
Pressure Detection:	Capacitive pressure sensor
Measurement Methods:	Oscillometric method
Pulse Wave Detection:	Capacitance type pressure sensor
Power Source:	4 "AA" batteries (not included) or AC adapter (sold separately Model HEM-ADPT1)
Battery Life:	Approx. 300 usages
Operating Temperature/Humidity:	50°F to 104°F (10°C to 40°C) 30 to 85% RH maximum
Storage Temperature/Humidity:	-4°F to 140°F (-20°C to 60°C) 10 to 95% RH maximum
Monitor Weight:	Approximately 13 oz (350 g) with cuff, without batteries
Outer Dimensions:	Approximately 5" (w) x 4 3/4" (d) x 3 3/4" (h) (125 mm x 119 mm x 95 mm)
Cuff Size:	Approximately 5 1/2" x 18 7/8" (140 mm x 480 mm)
Cuff Circumference:	Fits arm circumferences 9" to 13" (229 mm x 330 mm)
Advanced Diagnostics	Irregular Heartbeat Detector, 30 memory with date/time stamp, Advanced Averaging
Contents:	Monitor, Arm Cuff, Instruction Manual
Accessories (not included):	Large Cuff (13" - 17") Model H-003D and AC Adapter Model HEM-ADPT1

Specifications and appearance are subject to change without notice for purposes of improvement.

OMRON®
OMRON HEALTHCARE, INC.
1200 Lakeside Drive
Bannockburn, Illinois 60015
Telephone: 1-847-680-6200
Toll Free Order Number: 1-800-231-4030
Fax Number: 1-800-637-6763
Visit our Web site at www.omronhealthcare.com

OMRON HEALTHCARE CANADA
5500 N. Service Road Penthouse
Burlington, Ontario L7L 6W6
Telephone: 1-905-315-8408
Toll Free Order Number: 1-800-231-4030
Fax Number: 1-800-637-6763
Visit our Web site at www.omronhealthcare.com

© 11/2005 Omron Healthcare, Inc.

Printed in U.S.A.
712CN2PS Rev.