

**#1** Brand Recommended By Pharmacists and Doctors

**OMRON®**

HEM-712CLC


# Automatic Blood Pressure Monitor with Large Cuff



American Heart  
Association®   
*Learn and Live.™*



**New &  
Improved**

-  **IntelliSense®** Technology for More Accurate, More Comfortable Readings
- **Simple, Silent, 1-Touch Automatic Operation**
- **Extra Large Digital Display**
  - Systolic and Diastolic Blood Pressure
  - Pulse
  - Date and Time
- **Advanced Diagnostics**
  - Irregular Heartbeat Detector (♡)
  - 30 Memory Storage (📄)
  - Advanced Averaging Displays Up to Last 3 Readings in Last 10 Minutes (📊)
- **Durable Contoured Arm Cuff**
  - Large Cuff Fits Arms 13"-17" in Circumference
  - D-ring Cuff is Easy to Apply and Fits Securely
- **AC Adapter, Carrying Case, and Illustrated Instruction Manual Included (batteries not included)**
- **2 Year Warranty**

**A Good Sense of Health**®

# MODEL HEM-712CLC

## PACKAGING SPECIFICATIONS...

CARTON:	QTY:	APPROX. WEIGHT:	APPROX. DIMENSIONS: (W x D x H)
Unit: Package	1 each	2.25 lbs.	5.25" x 4.87" x 7"
Master:	12 each	30.2 lbs.	20.35" x 17.83" x 8.18"
UPC:	0 73796 07124 0		

## PRODUCT SPECIFICATIONS...

Display:	Liquid Crystal Display
Measurement Range:	Pressure: 0 - 280 mmHg Pulse: 40 - 180/minute
Accuracy/Calibration:	Pressure: $\pm 3$ mmHg or 2 % of reading Pulse: $\pm 5\%$ of reading
Inflation:	Automatic by electric pump
Deflation:	Automatic pressure release valve
Rapid Air Release:	Automatic exhaust valve
Pressure Detection:	Capacitive pressure sensor
Measurement Methods:	Oscillometric method
Pulse Wave Detection:	Capacitance type pressure sensor
Power Source:	4 "AA" batteries (not included) or AC adapter
Battery Life:	Approx. 300 usages
Operating Temperature/Humidity:	50°F to 104°F (10°C to 40°C) 30 to 85% RH maximum
Storage Temperature/Humidity:	-4°F to 140°F (-20°C to 60°C) 10 to 95% RH maximum
Monitor Weight:	Approximately 13 oz (350 g) without batteries
Outer Dimensions:	Approximately 5" (w) x 4 3/4" (d) x 3 3/4" (h) (125 mm x 119 mm x 94 mm)
Cuff Size:	Approximately 6 7/8" x 24 1/2" (175 mm x 622 mm)
Cuff Circumference:	Fits arm circumferences 13" to 17" (330 mm x 432 mm)
Advanced Diagnostics	Irregular Heartbeat Detector, 30 memory with date/time stamp, Advanced Averaging
Contents:	Monitor, Large Arm Cuff, Instruction Manual, AC Adapter, Carrying Case

Specifications and appearance are subject to change without notice for purposes of improvement.

**OMRON**®  
OMRON HEALTHCARE, INC.  
1200 Lakeside Drive  
Bannockburn, Illinois 60015  
Telephone: 1-847-680-6200  
Toll Free Order Number: 1-800-231-4030  
Fax Number: 1-800-637-6763  
Visit our Web site at [www.omronhealthcare.com](http://www.omronhealthcare.com)

OMRON HEALTHCARE CANADA  
5500 N. Service Road Penthouse  
Burlington, Ontario L7L 6W6  
Telephone: 1-905-315-8408  
Toll Free Order Number: 1-800-231-4030  
Fax Number: 1-800-637-6763  
Visit our Web site at [www.omronhealthcare.com](http://www.omronhealthcare.com)

© 11/2005 Omron Healthcare, Inc.

Printed in U.S.A.  
712CLCN2PS Rev.