

OVER-BED TABLES



		Length	Width
TS-160	top	15"	31 1/2"
	base	16"	26 1/2"
	height range	31" – 41 1/2"	

U-Base, Over Bed Table

- Walnut laminate top
- Chrome plated base
- U-base for easy access
- Dual wheel nylon casters
- Single hand, spring loaded height adjustments
- 50 lbs. load capacity
- Shipping weight: 29 lbs.

TS-165

U-Base, Over-Bed Table

- Same as TS-160 with Gray laminate top



		Length	Width
TS-199	top	15"	31 1/2"
	base	16"	26 1/2"
	height range	20"– 27 1/4"	

Low Height, H-Base, Over-Bed Table

- Walnut laminate top
- Chrome plated base
- H-Base design
- Dual wheel nylon casters
- Single hand, spring loaded height adjustments
- 50 lbs. load capacity
- Shipping weight: 29 lbs.



		Length	Width
TS-170	top	15"	31 1/2"
	base	16"	26 1/2"
	height range	31" – 41 1/2"	

H-Base, Over-Bed Table

- Walnut laminate top
- Chrome plated base
- Dual wheel nylon casters
- Single hand, spring loaded height adjustments
- 50 lbs. load capacity
- Shipping weight: 29 lbs.

TS-175

H-Base, Over-Bed Table

- Same as TS-170 with Gray laminate top



		Length	Width
TS-180	top	15"	31 1/2"
	base	16"	26 1/2"
	height range	32 1/2" – 43"	

H-Base, Tilt-Top, Over-Bed Table

- Walnut laminate top
- Chrome plated base
- Tilting top for easy patient use
- 1/4" lip edge on long sides
- Dual wheel nylon casters
- Single hand, spring loaded height adjustments
- 25 lbs. load capacity
- Shipping weight: 29 lbs.

TS-185

H-Base Tilt-Top, Over-Bed Table

- Same as TS-180 with Gray laminate top

		Length	Width
TS-190	top	15"	31 1/2"
	base	16"	26 1/2"
	height range	32 1/2" – 43"	

U-Base Tilt-Top Over-Bed Table

- Walnut laminate top
- Chrome plated base
- U-base for easy access
- Tilting top for easy patient use
- 1/4" lip edge on long sides
- Dual wheel nylon casters
- Single hand, spring loaded height adjustments
- 25 lbs. load capacity
- Shipping weight: 29 lbs.

TS-195

U-Base, Tilt-Top, Over-Bed Table

- Same as TS-190 with Gray laminate top