

TREATMENT CABINETS



	Length	Height	Depth
8920	22"	33"	17"

Mobile Equipment Cart with 2 Drawers

- Mounted on a stable 5-leg black nylon base
- 2 Euro-style drawers
- Top rimmed on 3 sides
- All laminate cabinet



	Length	Height	Depth
8940	22"	33"	17"

Mobile Equipment Cart with 1 Drawer and 1 Shelf

- Mounted on a stable 5-leg black nylon base
- 1 Euro-style drawer
- Top and back have wire feed cut-outs
- Top rimmed on 3 sides
- All laminate cabinet



	Length	Height	Depth
089	17 1/2"	5"	6 3/4"

Optional Basket

- Wire basket with white epoxy finish



	Length	Height	Depth
8990	44"	32"	24"

Mobile Splinting Cart

- Grey laminate construction
- 1 door, 1 drawer and 2 shelves storage
- Vertical dividers for casting material
- 2 push-handles
- 4" easy-roll, locking casters
- Top provides ample room for casting pan
- Extra tall countertop height

