

Pediatric Beam Scale

Mechanical pediatric scales provide precision accuracy and ensure child safety. Durable steel bases and easy-to-clean trays and seats guarantee long-lasting reliability.

Features include:

- Heavy-duty steel base for improved strength and durability
- Rotating poise bars for dual reading.

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

www.healthometermedical.com



Health o meter®
Professional

Pediatric Beam Scale

Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place

Bridgeview, Illinois 60455

USA Telephone: +1 (708) 598-9100

Facsimile: +1 (708) 233-5472

Customer Service

US and Canada

Telephone: (800) 815-6615

Facsimile: (800) 654-7330

HomProCS@pelstarllc.com

Customer Service

outside the US and Canada

Telephone: (708) 598-9100

Facsimile: (708) 599-0150

HomProCS@pelstarllc.com



- Heavy-duty steel base for improved strength and durability
- Sanitary tray with built-in measuring tape provides safety, reads up to 23" x 1/4" or 58 cm x 1 cm. (1522KL). Tray made from easy-to-clean ABS plastic
- Child seat provides safety and comfort and reduces infant movement (1524KL). Made from easy-to-clean ABS plastic
- Rotating poise bars for dual reading of pounds and kilograms. Poise bar lock-in for error-free operation
- Factory calibrated to ensure accuracy
- 2 Year Limited Warranty

Product Specifications

Model	Capacity	Increment	Product Wt.	Seat/Tray Dims	Product Dims
1522KL	130 lb/65 kg	1 oz/20 g	29 lb	14 1/4" x 24" x 2 1/2"	24" x 20 1/4" x 22 1/4"
1524KL	130 lb/65 kg	1 oz/20 g	36 lb	25" x 13" x 6 3/4"	25" x 20 1/4" x 22 1/4"

Packaging Specifications

Model	UPC	Case Pack	I2 of 5 Code	Master Carton Dims	Case Cube	Package Wt.
1522KL	027045642361	1	UPC	13" x 21" x 26"	4.11	37 lb
1524KL	027045642378	1	UPC	18 1/4" x 20 3/4" x 25 1/2"	5.59	45 lb

Health o meter®
Professional