

Pro Series Large Raised Dial, Large Platform Scales

With superior accuracy, this scale gives precise weight measurements. Extra large platform accomodates multiple users.

Features include:

- Heavy-duty steel base and platform for extra strength and durability
- Professional 8" easy-to-read dial
- Extra large platform for stability

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

www.healthometermedical.com



160LB
(Pounds only)



160KG
(Kilograms only)



160KL
(Pounds and kilograms)

Health o meter®
Professional

Pro Series Large Raised Dial, Large Platform Scales

Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place

Bridgeview, Illinois 60455

USA Telephone: +1 (708) 598-9100

Facsimile: +1 (708) 233-5472

Customer Service

US and Canada

Telephone: (800) 815-6615

Facsimile: (800) 654-7330

HomProCS@pelstarllc.com

Customer Service

outside the US and Canada

Telephone: (708) 598-9100

Facsimile: (708) 599-0150

HomProCS@pelstarllc.com



160LB
(Pounds only)



160KG
(Kilograms only)



160KL
(Pounds and kilograms)

- Raised dial for extra foot room
- Non-slip, textured mat
- Heavy-Duty steel base and platform for improved strength and durability
- Professional 8" easy-to-read dial
- Extra large platform for stability
- 1 Year Limited Warranty

Product Specifications

Model	Capacity	Increment	Product Wt.	Platform Dims	Dial Size
160LB	400 lb	1 lb	9 lb	18 1/2" x 11 1/2"	8" diameter
160KL	400 lb/180 kg	1 lb/1 kg	9 lb	18 1/2" x 11 1/2"	8" diameter
160KG	180 kg	1 kg	9 lb	47 cm x 29 cm	20 cm diameter

Packaging Specifications

Model	UPC	Case Pack	12 of 5 Code	Master Carton Dims	Case Cube	Package Wt.
160LB	080006108481	2	40080006108489	17 1/4" x 21 1/2" x 15"	3.22	24 lb
160KL	080006108474	2	40080006108472	17 1/4" x 21 1/2" x 15"	3.22	24 lb
160KG	080006108467	2	40080006108465	44 cm x 55 cm x 38 cm	3.22	24 lb

Health o meter®
Professional