

ProPlus® Digital Stretcher/Lift Scale

Stretcher scale with motion detection system to read actual weight. Industrial strength construction for long-lasting reliability.

- Hoyer lift
- 180° swivel head for side read
- Display pivots for easy reading by operators of all sizes

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

www.healthometermedical.com



Health o meter®
Professional

ProPlus® Stretcher/Lift Scale



Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place

Bridgeview, Illinois 60455

USA Telephone: +1 (708) 598-9100

Facsimile: +1 (708) 233-5472

Customer Service

US and Canada

Telephone: (800) 815-6615

Facsimile: (800) 654-7330

HomProCS@pelstarllc.com

Customer Service

outside the US and Canada

Telephone: (708) 598-9100

Facsimile: (708) 599-0150

HomProCS@pelstarllc.com



Also available

Disposable stretcher covers
(Model 63807)

2000KLS

- 270 patient memory
- 3-way tare with memory
- Updatable firmware
- Casters with wheel lock for mobility and stability.
- Customizable 1" backlit LCD displays in fractions and decimals
- Includes 1 USB and 2 serial ports
- Keypad entry, BMI, manual tare
- 6 D-cell batteries (not incl.) 120V AC adapter (incl.)
- Metallic gray

Product Specifications

Model	Capacity	Increment	Shipping Wt.	Stretcher Dims	Product Dims
2000KL	400 lb/180 kg	0.2 lb/1/4 lb/ 4 oz/0.1 kg	145 lb	61" (length)	45" x 26" x 72" (fully extended)

Packaging Specifications

Model	Case Pack	Master Carton Dims	Case Cube	Package Wt.
2000KL	1	2 cartons: 51 1/2" x 25 1/2" x 13 20" x 47" x 10"	9.88 5.44	79 lb 66 lb

Health o meter®
Professional