




Electronic Low Profile Veterinary Scale

Animals of any size can be weighed safely and accurately with Health o meter's veterinary scale. Auto-zero, auto-tare and weight lock-in offer supreme precision helping animal doctors provide the best care.

Features include:

- Stainless steel platform with heavy-duty, removable rubber mat
- Auto-zero, auto-tare, weight lock-in for readings after animal is removed, auto-reweigh (no need to remove animal from scale to reweigh)



2842KL

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

www.healthometermedical.com

Health o meter®
P r o f e s s i o n a l

Electronic Low Profile Veterinary Scale



Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place

Bridgeview, Illinois 60455

USA Telephone: +1 (708) 598-9100

Facsimile: +1 (708) 233-5472

Customer Service

US and Canada

Telephone: (800) 815-6615

Facsimile: (800) 654-7330

HomProCS@pelstarllc.com


Customer Service

outside the US and Canada

Telephone: (708) 598-9100

Facsimile: (708) 599-0150

HomProCS@pelstarllc.com



- Motion-sensing weighing technology
- 1" Remote LCD display with bracket for wall-mount or table-top display
- Heavy duty steel base
- Low profile (2") 42" x 22 3/8" oversize platform
- Removable rubber mat for easy cleaning
- Reweigh, LB/KG, On/Off and Zero buttons
- Low battery indicator, auto-off feature
- 6 C-cell batteries & 120V AC adapter
- 1 Year Warranty

Product Specifications

Model	Capacity	Increment	Product Wt.	Platform Dims	Display Dims
2842KL	600 lb / 270 kg	0.2 lb/0.1 kg	40 lbs	42" x 22 3/8"	7" x 4 1/4" x 7 1/4"

Packaging Specifications

Model	UPC	Case Pk	I2 of 5 Code	Master Carton Dims	Case Cube	Pkg Wt.
2842KL	027045642491	1	UPC	50" x 24" x 5.5"	3.82	45 lb

Health o meter®
Professional