



Remote Display Scales

These high precision floor scales provide the added benefit of remote displays that can be mounted on a wall or placed on a table for easy reading. Steel construction with a large stable platform guarantees safety and durability.

Features include:

- Remote LCD display on 9-foot cord
- Weighs in kilograms and pounds
- 18 gauge steel construction with baked enamel finish

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

www.healthometermedical.com



320KL



349KLX

Health o meter®
Professional

Remote Display Scales



Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place
 Bridgeview, Illinois 60455
 USA Telephone: +1 (708) 598-9100
 Facsimile: +1 (708) 233-5472

Customer Service

US and Canada

Telephone: (800) 815-6615
 Facsimile: (800) 654-7330
HomProCS@pelstarllc.com

Customer Service

outside the US and Canada

Telephone: (708) 598-9100
 Facsimile: (708) 599-0150
HomProCS@pelstarllc.com



320KL



349KLX

- 0.2 lb/0.1 kg graduation is the best available in a floor scale (349KLX)
- Remote LCD display on 9' cord may be mounted on wall or placed on a table using display bracket (included) for ease of reading
- Weight is displayed in pounds or in kilograms
- 18 gauge steel construction with baked enamel finish
- Tare feature and 30-second Auto-Off to conserve energy. Low-battery indicator.
- 349KLX includes RS232 serial port for computer interface
- 9 volt battery and AC adapter included
- Non skid feet on a large 12" x 12" platform with non-skid mat
- 320KL – silver, 349KLX – white
- 1 Year Limited Warranty

Product Specifications

Model	Capacity	Increment	Product Wt.	Platform Dims	Product Dims
320KL	400 lb/180 kg	0.5 lb/0.2 kg	9 lb	12" x 12"	12" x 12" x 21/2"
349KLX	400 lb/180 kg	0.2 lb/0.1 kg	9 lb	12" x 12"	12" x 12" x 21/2"

Packaging Specifications

Model	UPC	Case Pk	12 of 5 Code	Master Carton Dims	Case Cube	Pkg Wt.
320KL	080006108535	2	40080006108533	6 1/2" x 19" x 13 1/2"	0.96	24 lb
349KLX	080006108542	2	40080006108540	6 1/2" x 19" x 13 1/2"	0.96	24 lb

Health o meter®
 Professional