

Portable Veterinary Mechanical Scale with Built-In Tray

Our portable veterinary scales let you take strength and accuracy from room to room. The dial is angled for easy reading and the sanitary tray provides extra safety. Just grab the convenient carrying handle and you're on your way.

Features include:

- Convenient carry handle for easy transport
- Easy-to-read angled dial
- Weight lock for fool-proof reading after weight is removed

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

www.healthometermedical.com



386S-01

Health o meter®
Professional

Portable Veterinary Mechanical Scales



Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place

Bridgeview, Illinois 60455

USA Telephone: +1 (708) 598-9100

Facsimile: +1 (708) 233-5472

Customer Service

US and Canada

Telephone: (800) 815-6615

Facsimile: (800) 654-7330

HomProCS@pelstarllc.com

Customer Service

outside the US and Canada

Telephone: (708) 598-9100

Facsimile: (708) 599-0150

HomProCS@pelstarllc.com



386S-01

- Portable veterinary scale with convenient carrying handle for easy transport
- Easy-to-read angled dial
- Weight Lock for error-proof reading after weight is removed.
- Sanitary tray provides safety. Tray is removable, made from easy-to-clean ABS plastic.
- Measuring tape (18"/46 cm) included on tray
- Protected zero adjust knob for easy access and to eliminate breakage.
- 1 Year Limited Warranty

Product Specifications

Model	Capacity	Increment	Product Wt.	Tray Dims	Product Dims
386S-01	50 lb	1/4 lb	10 lb	19 1/4" x 12 3/4" x 2"	19 3/8" x 12 1/4" x 5 1/2"
386KGS-01	25 kg	100 g	10 lb	19 1/4" x 12 3/4" x 2"	19 3/8" x 12 1/4" x 5 1/2"

Packaging Specifications

Model	UPC	Case Pack	I2 of 5 Code	Master Carton Dims	Case Cube	Package Wt.
386S-01	080006009801	1	20080006009805	21 1/2" x 6 1/2" x 16 1/4"	1.31	12 lb
386KGS-01	080006009818	1	20080006009812	21 1/2" x 6 1/2" x 16 1/4"	1.31	12 lb

Health o meter®
Professional