

Portable Dial Seat Scale

When patients are unsteady on their feet, these scales help you take the weight measurements you need. Their comfortable seats and ability to lock in weights make them perfect for home care, nursing home or hospital use.

Features include:

- Easy-to-read angled dial
- Weight lock for fool-proof reading after weight is removed
- Protected zero adjust knob for easy access and to eliminate breakage

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

www.healthometermedical.com



388S-01

Health o meter®
Professional

Portable Seat Scale



Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place

Bridgeview, Illinois 60455

USA Telephone: +1 (708) 598-9100

Facsimile: +1 (708) 233-5472

Customer Service

US and Canada

Telephone: (800) 815-6615

Facsimile: (800) 654-7330

HomProCS@pelstarllc.com

Customer Service

outside the US and Canada

Telephone: (708) 598-9100

Facsimile: (708) 599-0150

HomProCS@pelstarllc.com



388S-01

- Portable seat scale perfect for home care, nursing home and hospital use
- Easy-to-read angled dial
- Weight Lock for error-proof reading after patient is removed
- Protected zero adjust knob for easy access and to eliminate breakage.
- Contoured seat for balance and comfort
- Convenient carrying handle for portability
- 1 Year Limited Warranty

Product Specifications

Model	Capacity	Increment	Product Wt.	Seat Dims	Product Dims
388S-01	250 lb	1 lb	9 lb	14 1/4" x 11 1/4"	15 1/2" x 12 1/4" x 5 3/8"

Packaging Specifications

Model	UPC	Case Pack	I2 of 5 Code	Master Carton Dims	Case Cube	Package Wt.
388S-01	080006004110	1	20080006004114	18 1/4"W x 7"H x 15"D	1.11	13 lb

Health o meter®
Professional