

## Physician Balance Beam Scales

Trusted across the nation, our beam scales offer the industry standard in strength and durability. Doctors appreciate our accuracy and long-lasting reliability.

### Features include:

- Rotating poise bars for dual reading in pounds and kilograms. Poise bar lock-in for error free operation
- Permanent, engraved readings on poise bar for long lasting readability

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

[www.healthometermedical.com](http://www.healthometermedical.com)



Health o meter®  
Professional

# Physician Beam Scales



Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place

Bridgeview, Illinois 60455

USA Telephone: +1 (708) 598-9100

Facsimile: +1 (708) 233-5472

Customer Service

US and Canada

Telephone: (800) 815-6615

Facsimile: (800) 654-7330

[HomProCS@pelstarllc.com](mailto:HomProCS@pelstarllc.com)

Customer Service

outside the US and Canada

Telephone: (708) 598-9100

Facsimile: (708) 599-0150

[HomProCS@pelstarllc.com](mailto:HomProCS@pelstarllc.com)



402KL



402LB



400KL

- Rotating poise bars for dual reading of pounds and kilograms (402KL, 400KL). Poise bar lock-in for error-free operation
- Permanent, engraved readings on poise bar for long-lasting readability
- Heavy-duty steel base for improved strength and durability
- Professional “Silent Slide™” height rod for smooth, quiet operation. Dual-read height rod (402KL, 402LB) range 24" to 84" (by 1/8 in.), 60cm to 213cm (by 1 mm).
- ABS plastic platform cover for easy cleaning
- Factory calibrated to ensure accuracy
- Handbook 44 certifiable (402LB only)
- 402LB – Includes 100 lb counterweight to increase capacity to 450 lbs
- 402KL – Includes 100 lb and 40 kg counterweights to increase capacity to 450 lbs/ 200 kgs
- Available Accessories: Counterweights increase capacity to 450 lbs. or 200kg, casters provide easy mobility
- Assembly required
- 2 Year Limited Warranty

## Product Specifications

Model	Capacity	Increment	Product Wt.	Platform Dims	Product Dims
402KL	500 lb/200 kg	1/4 lb/100 g	39 lb	10 5/8" x 14"	10 5/8" x 20 1/2" x 57 1/4"
402LB	500 lb	1/4 lb	38 lb	10 5/8" x 14"	10 5/8" x 20 1/2" x 57 1/4"
400KL	350 lb/160 kg	1/4 lb/100 g	33 lb	10 5/8" x 14"	10 5/8" x 20 1/2" x 57 1/4"
402KLWH	500 lb/200 kg	1/4 lb/100 g	33 lb	10 5/8" x 14"	10 5/8" x 20 1/2" x 57 1/4"
402KLBWH	500 lb	1/4 lb/100 g	33 lb	10 5/8" x 14"	10 5/8" x 20 1/2" x 57 1/4"
402KLBCERT	500 lb	1/4 lb	33 lb	10 5/8" x 14"	10 5/8" x 20 1/2" x 57 1/4"

## Packaging Specifications

Model	UPC	Case Pack	I2 of 5 Code	Master Carton Dims	Case Cube	Package Wt.
402KL	027045642354	1	UPC	20 1/2" x 6 3/4" x 59 1/4"	4.74	45 lb
402LB	027045642408	1	UPC	20 1/2" x 6 3/4" x 59 1/4"	4.74	45 lb
400KL	027045642347	1	UPC	20 1/2" x 6 3/4" x 59 1/4"	4.74	40 lb
402KLWH	tbd	1	UPC	20 1/2" x 6 3/4" x 59 1/4"	4.74	45 lb
402KLBWH	tbd	1	UPC	20 1/2" x 6 3/4" x 59 1/4"	4.74	45 lb
402KLBCERT	022319408292	1	UPC	20 1/2" x 6 3/4" x 59 1/4"	4.74	40 lb

Health o meter®  
Professional