

Digital Floor Scales

These high precision floor scales give you accurate weight readings at home. They are easy to clean and easy to read with measurements in pounds or kilograms.

Features include:

- Easy-to-read super bright display
- Weight in kilograms and pounds
- All metal base with easy-to-clean plastic platform

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

www.healthometermedical.com



830KL



845KL



849KL

Health o meter®
Professional

Digital Floor Scales

Health o meter®
Professional Products
Pelstar LLC
Corporate Headquarters
 7400 West 100th Place
 Bridgeview, Illinois 60455
 USA Telephone: +1 (708) 598-9100
 Facsimile: +1 (708) 233-5472

Customer Service
US and Canada
 Telephone: (800) 815-6615
 Facsimile: (800) 654-7330
HomProCS@pelstarllc.com

Customer Service
outside the US and Canada
 Telephone: (708) 598-9100
 Facsimile: (708) 599-0150
HomProCS@pelstarllc.com



.8" LED

830KL



1.5" LED

845KL



1.5" LCD

849KL

- Easy to read displays
- Weigh in kilograms and pounds
- All metal base with easy to clean plastic platform
- Large stable platform
- 9V battery included (830KL, 845KL)
- Long life replaceable lithium battery (849KL only)
- **1 Year Limited Warranty**

Product Specifications

Model	Capacity	Increment	Product Wt.	Platform Dims	Product Dims
830KL	300 lb/136 kg	0.5 lb /0.2 kg	4 lb	12" x 10 1/2"	12 1/4" x 1 1/4" x 11"
845KL	330 lb/150 kg	0.2 lb/0.1 kg	7 lb	13 1/2" x 11 1/2"	13 1/2" x 2" x 11 1/2"
849KL	330 lb/150 kg	0.2 lb/0.1 kg	7 lb	13 1/2" x 11 1/2"	13 1/2" x 2" x 11 1/2"

Packaging Specifications

Model	UPC	Case Pack	I2 of 5 Code	Master Carton Dims	Case Cube	Package Wt.
830KL	080006108504	2	40080006108502	5 1/8" x 13 1/2" x 13 3/4"	0.55	10 lb
845KL	080006108511	2	40080006108519	6 1/4" x 15 1/4" x 13 1/2"	0.74	16 lb
849KL	080006108528	2	40080006108526	6 1/4" x 15 1/4" x 13 1/2"	0.74	16 lb

Health o meter®
Professional