

# MC-301N Basal Digital Thermometer



- Safe and accurate
- Helps identify ovulation timing
- Beeps at peak temperature
- Last temperature recall
- Includes temperature tracking chart

Lifetime Warranty

# MODEL MC-301N

## PACKAGING SPECIFICATIONS...

CARTON:	QTY:	APPROX. WEIGHT:	APPROX. DIMENSIONS: (l x w x h)	CASE CUBE
Unit:	1 Each	1.6 oz.	6.5" x 1.0" x 2.63"	.01
Master:	48 Each	6 lbs. 3 oz.	13 1/4" x 8 3/4" x 12"	.81

## PRODUCT SPECIFICATIONS...

Display:	LCD
Measurement Range:	90°F to 109.9°F 32.2°C to 43.2°C
Accuracy:	+/- 0.2°F between 96.0°F to 106.0°F (+/- 0.1°C between 35.6°C to 41.1°C). +/- 0.4°F at less than 96.0°F and greater than 106.0 °F (+/- 0.2°C at less than 35.6°C and greater than 41.1°C).
Response Time:	60 seconds
Display Resolution:	0.1°F
Power Source (one battery- 1.5v button type):	SR41
Battery Life (continuous use):	350 hrs.
Automatic Shut-Off (approximately):	10 min.
Operating Temperature/Humidity:	50°F to 104°F (10°C to 40°C) 15 to 95% R.H. noncondensation
Storage Temperature/Humidity:	-4°F to 140°F (-20C to 60°C)
Thermometer Weight with Battery:	.4 oz.
Thermometer Dimension (l x w x h):	5" x 7/16" x 3/8"
Waterproof:	Yes
UPC Code:	0 73796 03301 9

Specifications and appearance are subject to change without notice for purposes of improvement.