



Quality Medical Products Since 1919

Puritan Medical Products Company LLC  
Guilford, Maine 04443-0149 USA

An affiliate of Hardwood Products Company LLC

Tel 800-321-2313 (US and Canada)  
207-876-3311

Fax 800-323-4153 (US and Canada)  
207-876-3130

Web [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

Email [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

**MEDICAL WOODENWARE** Our wood-based products are made from Northern White Birch, which yields a high tensile strength while being smooth to the touch. Medical woodenware is used for patient oral exams, gathering specimens, applying and mixing medicine, and various health and beauty applications.

## APPLICATOR STICKS

Choose the right size and shape for your applications. Utilize for specimen transfer, as mixing rods, for laboratory applications, and other general-purpose functions. Made with straight-cut ends.

	Item No.	Description	Length & Dia.	Packing	Weight	
	Sterile 25-807 2W	Wood, sterile	6" x 1/12"	10/100/2 per pkg.	6 lbs.	
		805	Wood	12" x 1/12"	12/1000	21 lbs.
	Non-Sterile 807	Wood	5-3/4" x 1/12"	20/1000	18 lbs.	
		809	Wood	12" x 1/8"	12/500	23 lbs.
		811	Wood	18" x 1/8"	2/1000	12 lbs.

Products are not shown actual size.