



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
 Guilford, Maine 04443-0149 USA

An affiliate of Hardwood Products Company LLC

Tel 800-321-2313 (US and Canada)
 207-876-3311

Fax 800-323-4153 (US and Canada)
 207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

CYTOLOGY COLLECTION DEVICES

Puritan cytology products are made to maximize cell collection, while enhancing patient comfort. The Puritan CerviSoft® offers gentle specimen sampling for Pap smear tests.

	Item No.	Description	Shaft Length	Packing	Weight	
	25-2188	Histobrush®, sterile	8"	10/100/1 per pkg.	7 lbs.	
		25-2197	CerviSoft® Foam, sterile	7"	10/100/1 per pkg.	7 lbs.
		25-2199	Histobrush®, sterile	7"	10/100/1 per pkg.	7 lbs.
      	2188	Histobrush®	8"	10/100	6 lbs.	
	2194	Rovers® Spatula	7-3/4"	15/50	4 lbs.	
	2195	Rovers® Cervex-Brush®	8"	35/25	5 lbs.	
	2195 Combi	Rovers® Cervex-Brush® Combi	8"	4/25	2 lbs.	
	2196	Rovers® EndoCervex-Brush®	8"	10/50	4 lbs.	
	2197	CerviSoft® Foam	7"	10/100	2.5 lbs.	
	2199	Histobrush®	7"	10/100	6 lbs.	

Products are not shown actual size.

ENDOCERVEX-BRUSH® | 2196

1 Detachable tip allows use with both liquid-based and conventional cervical cytology tests

2 Patient-friendly, one-piece brush provides greater comfort than traditional wire and nylon-bristle design

3 Special safety ring decreases the chance of the brush entering the endometrium during an examination