



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

An affiliate of Hardwood Products Company LLC

Tel 800-321-2313 (US and Canada)
207-876-3311







Fax 800-323-4153 (US and Canada)
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

FOAM-TIPPED APPLICATORS

These applicators are ideal for a variety of clinical and DNA applications. Tips are made of medical-grade foam. Choose from a range of shapes and shaft materials.

	Item No.	Description	Shaft Length & Dia.	Packing	Weight
	25-1406 1PF	Ultrafine polystyrene shaft, sterile	6" x .095"/.052"	1000 bulk	5 lbs.
	25-1506 1PF	Plastic shaft, sterile	6" x 1/10"	10/100/1 per pkg.	5 lbs.
	25-1506 1PF TT	Plastic shaft, sterile, dry transport tube	6" x 1/10"	24/30/1 per tube	5 lbs.
	25-1605 1PSF RECT	Polypropylene shaft, rectangular tip, sterile	5" x .906"	10/50/1 per pkg.	4 lbs.
	25-1803 1PF	Polypropylene shaft, sterile	3" x 3/32"	10/100/1 per pkg.	5 lbs.
	25-1805 1PF RND	Polypropylene shaft, round tip, sterile	5" x .906"	10/50/1 per pkg.	6 lbs.

Products are not shown actual size.