



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

An affiliate of Hardwood Products Company LLC

Tel 800-321-2313 (US and Canada)
207-876-3311

Fax 800-323-4153 (US and Canada)
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

POLYESTER-TIPPED APPLICATORS

Synthetic applicator tips are made of medical-grade polyester that we card ourselves. They're non-absorbent and contain no fatty acids or latent sugars, making them ideal in laboratory settings and for inclusion in many diagnostic kits.

	Item No.	Description	Shaft Length & Dia.	Packing	Weight
Sterile	25-800 D	Aluminum shaft, sterile	5-1/2" x .035"	10/100/1 per pkg.	5 lbs.
	25-800 1PD	Ultrafine plastic shaft, sterile	6" x .095"/.052"	10/100/1 per pkg.	5 lbs.
	25-800 1PD ALUM	Plastic/aluminum shaft, sterile	6" x .035"	10/100/1 per pkg.	5 lbs.
	25-801 D	Aluminum shaft, sterile	5-1/2" x .025"	10/100/1 per pkg.	5 lbs.
	25-806 1PD	Plastic shaft, sterile	6" x 1/10"	10/100/1 per pkg.	5 lbs.
	25-806 2PD			10/100/2 per pkg.	6 lbs.
	25-806 1WD	Wood shaft, sterile	6" x 1/10"	10/100/1 per pkg.	5 lbs.
	25-806 2WD			10/100/2 per pkg.	6 lbs.
	25-806 1PP TT Cap	Plastic shaft, sterile, dry transport tube, molded end cap	6" x 1/10"	24/30/1 per tube	7 lbs.
	25-806 1PD TT	Plastic shaft, sterile, dry transport tube	6" x 1/10"	24/30/1 per tube	7 lbs.
25-806 2PD TT	24/30/2 per tube			7 lbs.	
25-1000 1PD	Plastic tapered shaft, sterile, scored	6" x .085"/.040"	10/100/1 per pkg.	5 lbs.	
Non-Sterile	803-PD	Plastic shaft	3" x 1/10"	10/10/100	9 lbs.
	806-PD	Plastic shaft	6" x 1/10"	10/10/100	15 lbs.
	806-WD	Wood shaft	6" x 1/10"	10/10/100	12 lbs.
	1000-PD	Plastic tapered shaft, scored	6" x .085"/.040"	10/100	0.6 lbs.

Products are not shown actual size.

