



Quality Medical Products Since 1919

Puritan Medical Products Company LLC  
Guilford, Maine 04443-0149 USA

An affiliate of Hardwood Products Company LLC

Tel 800-321-2313 (US and Canada)  
207-876-3311

Fax 800-323-4153 (US and Canada)  
207-876-3130

Web [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

Email [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

## RAYON-TIPPED APPLICATORS

Ideal for diagnostic specimen collection, wound treatment, and application of medicines. Tips are made from the finest USP grades of rayon fiber, which we card in our Guilford facility.

	Item No.	Description	Shaft Length & Dia.	Packing	Weight
	25-800 R	Aluminum shaft, sterile	5-1/2" x .035"	10/100/1 per pkg.	5 lbs.
	25-801 R	Aluminum shaft, sterile	5-1/2" x .025"	10/100/1 per pkg.	5 lbs.
	25-806 1PAR	Fine paper shaft, sterile	6" x .059"	10/100/1 per pkg.	5 lbs.
	25-806 1PR	Plastic shaft, sterile	6" x 1/10"	10/100/1 per pkg.	5 lbs.
	25-806 1WR	Wood shaft, sterile	6" x 1/10"	10/100/1 per pkg.	5 lbs.
	25-808 2PR	Polypropylene shaft, throat or OB/GYN swab, sterile	8" x 5/32"	10/50/2 per pkg.	10 lbs.
	803-PR	Plastic shaft	3" x 1/10"	10/10/100	9 lbs.
	806-PR	Plastic shaft	6" x 1/10"	10/10/100	15 lbs.
	806-WR	Wood shaft	6" x 1/10"	10/10/100	12 lbs.
	808	Paper shaft, throat or OB/GYN swab	8" x 5/32"	10/50	5 lbs.
	816	Polypropylene shaft, proctoscopic applicator	16" x 5/32"	10/50	9 lbs.

Products are not shown actual size.