

ZOLL AED Plus Incident Log

34018 2002-05-22 15:43:06 X01K000217 348939 DOE, JANE R.

Patient

Patient ID / MR Number: 348939
Last name: DOE
First name: JANE
MI: R
Gen:
Sex: Female
Race:
Date of birth: 01/15/1904
Height:
Weight:

Case

Date: 05/22/2002
Start Time: 15:43:06
Run Number / CPR Number: 34018
Device Type: AED Plus
Serial number: X01K000217
Device ID:
Software Version: 1.11
Operator: Medic04
Case Length: 00:04:36
Number of Shocks Delivered: 5

Event Summary

Time:	Event:
15:43:43	Analysis: Started
15:43:52	Analysis: Shock Advised
15:43:56	Shock:1 120 Joules
15:43:56	Patient Impedance:56 ohms
15:44:01	Analysis: Started
15:44:10	Analysis: Shock Advised
15:44:17	Shock:2 150 Joules
15:44:17	Patient Impedance:56 ohms
15:44:22	Analysis: Started
15:44:31	Analysis: Shock Advised
15:44:37	Shock:3 200 Joules
15:44:37	Patient Impedance:55 ohms
15:45:50	Analysis: Started
15:45:59	Analysis: Shock Advised
15:46:05	Shock:4 200 Joules
15:46:05	Patient Impedance:56 ohms
15:46:10	Analysis: Started
15:46:19	Analysis: Shock Advised
15:46:25	Shock:5 200 Joules
15:46:25	Patient Impedance:55 ohms
15:46:30	Analysis: Started

ZOLL AED Plus Incident Log

34018 2002-05-22 15:43:06 X01K000217 348939 DOE, JANE R.

15:46:39

Analysis: No Shock Advised

Comments:

AED Plus sample data

Physician Signature _____

Print Name _____