



Welch Allyn
Green Series™ 777
Integrated Wall System



WelchAllyn®

Welch Allyn Green Series 777 Integrated Wall System

Allows for reliable diagnostic patient data in a standardized and energy efficient solution to improve your bottom line.

Welch Allyn

Green Series 777 Wall Transformer—configured to each facility's needs.

- Help reduce facility cost with a more energy efficient product than the 767 Wall Transformer
- Help improve workflow and save energy with handles that automatically turn on/off
- Quickly and easily upgrade from the Welch Allyn 767 Wall Transformer to the Green Series 777 Wall Transformer
- No on/off switch reduces contact, to help reduce the risk of cross-contamination



Welch Allyn Oscopes & Ophthalmoscopes

- Standardize or choose the technology you need for each facility
- Upgrade to the Welch Allyn LED lamp for greater energy efficiency and lower cost of ownership



Welch Allyn MacroView™ Otoscope

- Almost twice the field of view and 30% greater magnification than traditional otoscopes
- Allows for simplified diagnosis of otitis media through superior visualization of the tympanic membrane



Welch Allyn Pneumatic Otoscope

- Rotating lens with ultraseal enables easy pneumatic otoscopy
- Easier instrumentation to allow a faster and more comfortable exam for the patient



Welch Allyn PanOptic™ Ophthalmoscope

- Fast, easy entry into undilated pupils with 5X larger view of the fundus for a streamlined exam
- 25° field of view vs. the standard 5° field of view provides optimal visualization



Welch Allyn Coaxial Ophthalmoscope

- 18 unique aperture/filter combinations for greater versatility
- 28 focusing lenses with a range of -25 to +40 diopters



Welch Allyn industry-leading blood pressure and thermometry devices provide you with fast, accurate and reliable patient data you can trust—to help improve patient care and clinical decisions.



Welch Allyn Wall Aneroid

- Recessed, laser-engraved dial face with high-contrast pointer provides certified accuracy to reduce errors and improve patient care
- Lifetime calibration warranty reduces cost of ownership



Welch Allyn Connex® ProBP™ Digital Blood Pressure Device

- Connectivity to most leading EHR systems, to improve patient documentation accuracy and efficiency
- Improve patient comfort and increase throughput with fast, 15-second NIBP readings
- Portable, so you can take BP readings just about anywhere



Welch Allyn Spot Vital Signs® Device

- Improve workflow by using a device that includes NIBP, pulse rate and MAP
- Optional SureTemp® thermometry and Masimo® or Nellcor® pulse oximetry



Welch Allyn Connex® Spot Monitor

- Improve patient comfort and increase throughput with 15-second NIBP readings taken as the cuff inflates
- Automated blood pressure averaging helps improve hypertension diagnosis accuracy
- Connectivity to most leading EHR systems, to improve patient vitals documentation accuracy and efficiency
- Improve workflow by using a device that includes NIBP, pulse rate, MAP, thermometry and SpO₂



Welch Allyn SureTemp® Plus Electronic Thermometer

- Provides accurate oral temperatures in 4 to 6 seconds, axillary and rectal temperatures in 10 to 13 seconds and adult axillary temperatures in 12 to 15 seconds
- Reduce risk of cross-contamination with single-use probe covers
- Improve efficiency with convenient storage for 75 probe covers
- Takes approximately 6,000 readings on 3 AA batteries and provides last temperature recall

Standardize across facilities for greater efficiency and consistent data capture or *customize* for your facility's needs.

	WALL BOARD	GS 777 WALL TRANSFORMER	COAXIAL OPHTHALMOSCOPES	PANOPTIC™ OPHTHALMOSCOPES		DIAGNOSTIC MACROVIEW™ OTOSCOPIES		PNEUMATIC OTOSCOPE	KLEENSPEC® DISPENSERS		WALL ANEROID	CONNEX® PROBP™ DEVICE	SURETEMP® PLUS THERMOMETER
Part #		77710	11720	11810	11820	23810	23820	20200	52400-PF	52100-PF	7670-01	34XFWT-B	01690-300
GS 777 Integrated Wall Systems													
77798	•	•	•					•		•	•		
77791-1MPX	•	•		•		•			•		•		•
77791-2MPX	•	•			•	•			•		•		•
77791-MX	•	•	•			•			•		•		•
77792-1MP	•	•		•		•			•		•		
77792-2MP	•	•			•	•			•		•		
77792-M	•	•	•			•			•		•		
77795-2MPX	•	•			•		•		•		•		•
77795-MX	•	•	•				•		•		•		•
77796-M	•	•	•				•		•		•		
77797-2PX	•	•			•			•		•	•		•
77797-X	•	•	•					•		•	•		•
77791-MXNOBP	•	•	•			•			•			•	•
77792-MNOBP	•	•	•			•			•				•
GS 777 Integrated Wall Systems with Connex® ProBP™ Digital Blood Pressure Device (includes wall mounting bracket and device)													
77791-MXPROBP	•	•	•			•			•			•	•
77791-2MPXPROBP	•	•			•	•			•			•	•
GS 777 Integrated Wall Systems with Connex® Spot Monitor, Spot Vital Signs® or Spot Vital Signs® LXi Devices (includes wall mounting bracket only, devices must be ordered separately)													
77794-MCSM	•	•	•			•			•				
77794-2MCSM	•	•			•	•			•				
77794-MSPOT	•	•	•			•			•				
77794-MLXI	•	•	•			•			•				
77794-2MPSPOT	•	•			•	•			•				
77794-2MPLXI	•	•			•	•			•				
GS 777 Wall Transformer Sets													
77510		•							•				
77910		•							•		•		
77710-71M		•	•			•							
77710-81M		•		•		•							

The United States Green Building Council's Leadership in Energy and Environmental Design (LEED) is an internationally recognized green building certification system providing benchmarks for the design, construction and operation of high-performance green buildings.

Welch Allyn is proud to have earned a Gold rating in the LEED certification program.

To place an order, call our Customer Service Department at: 1.800.535.6663 or contact your local Welch Allyn representative.

© 2016 Welch Allyn MC13012 80021007 Ver A 2016-01 www.welchallyn.com

