

# FF28LWH



24.88" x 18.5" x 17.75" (H x W x D)

Compact all-refrigerator with automatic defrost, front-mounted lock, and white finish; Replaces FF28L

### Highlights:

Compact size is ideal for use on top of counters

Save time and energy with low maintenance automatic defrost

Factory-installed lock provides security you can count on

### Product Features:

Automatic defrost	Reduced maintenance with auto defrost system
Compact dimensions	Ideal size for countertop use
Fully finished cabinet	Allows the refrigerator to be used as a freestanding model
Factory installed lock	Keyed lock for secure interior
Flat door liner	Interior design conforms to CDC guidelines to ensure stable conditions for sensitive materials
Adjustable shelf	Rearrange refrigerator space to accommodate all sizes or remove shelf for simple cleanup
One piece interior liner	Easy clean-up with a one piece liner that won't hold a mess in hidden crevices
Reversible door	User reversible door swing
Adjustable thermostat	Manage the temperature with ease and accuracy

## FF28LWH Specifications:

Overview	
Height of Cabinet	24.5" (62 cm)
Height to Hinge Cap	24.88" (63 cm)
Width	18.5" (47 cm)
Depth	17.75" (45 cm)
Depth with door at 90°	34.5" (88 cm)
Capacity	2.4 cu.ft. (68 L)
Defrost Type	Automatic
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	1.3
Voltage/Frequency	115 V AC/60 Hz
Weight	36.0 lbs. (16 kg)
Shipping Weight	38.0 lbs. (17 kg)
Estimated Time to Ship	Next day shipping
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
Refrigerator Features	
Door Swing	RHD
Reversible	Yes
Shelf Type	Wire
Shelf Qty	1
Adjustable Shelves	Yes
Thermostat Type	Dial
Refrigerant Type	R134a
Refrigerant Amount	1.28oz.
High Side PSI	261.0
Low Side PSI	103.0
Level Legs Qty	2
Temperature Range	36°F to 46°F
Dimensions	
Interior Height	21.5" (55 cm)
Interior Width	16.0" (41 cm)
Interior Depth	13.5" (34 cm)
Compressor Step Height	7.75" (20 cm)
Compressor Step Width	16.0" (41 cm)

Compressor Step Depth 5.0"

