

SS8381

SS8300 Series Foot Stools



Stainless Steel Transport Cart



Shown with optional (UMF Medical item SS8380) Stackable

Standard Features

- High quality, welded stainless steel construction
- 16 gauge, 1" (2.54cm) outside diameter stainless steel tubing
- 1.25" x .125" (3.2cm x .3cm) stainless steel band below step platform for rigidity
- 2" (5.1cm) swivel casters
- Convenient push handle
- 3 year limited warranty

Specifications

- **Width:** 20.5" (52.1cm)
- **Height (platform):** 8.5" (21.6cm)
- **Height (push handle):** 36" (91.4cm)
- **Depth:** 15" (38.1cm)

