

## Flat Top, Classic Series, Straight Line Treatment Table with Shelf



### MODEL 200

- Solid hardwood legs
- H-brace construction
- Full laminate shelf in Natural or Dark Cherry to coordinate with frame
- Legs have pre-inserted carriage bolts for fast easy assembly
- Seamless, rounded corner top
- 2" of firm foam padding (5 cm)
- Available in Natural and Dark Cherry finish
- Choice of upholstery colors
- Fast, easy assembly with basic tools
- 400 lbs. load capacity under normal use (181.43 kg)



### SPECIFICATIONS

200-27	72" L x 27" W x 31" H
200-30	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:

**!** **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2

### Options



### Finishes



Natural

Dark Cherry

### Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black