

Alpha Series, Treatment Table with Shelving



MODEL 3030

- Features powder-coated, all-steel frame, gray laminate shelving and an adjustable backrest
- Easy bolt-through assembly
- Paper dispenser included
- Adjustable leg levelers
- Seamless rounded corner top
- 2" firm foam padding (5 cm)
- Heavyweight, knit-backed upholstery
- 450 lbs. load capacity under normal use (204 kg)

SPECIFICATIONS

3030-27	72" L x 27" W x 30" H
3030-30	72" L x 30" W x 30" H



California residents please be advised of the following, pursuant to Proposition 65:

! **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Options




020 Small Pillow



060 Large Pillow



009 Nosehole



077 Safety Strap



040 Paper Cutter

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black