



Hillrom™

A close-up photograph of a female doctor with dark, wavy hair, wearing a white lab coat and a pink stethoscope. She is holding a handheld medical device with a light, which she is using to examine the ear of a male patient whose face is partially visible on the left. The background is a blurred clinical setting.

**LET US
CHANGE
YOUR VIEW**

Explore Welch Allyn® Physical Exam Sets and Systems
featuring the new, LED wide-view PanOptic™ Plus
Ophthalmoscopes and MacroView® Plus Otoscopes

WALL SYSTEMS AND TRANSFORMERS

	OPHTHALMOSCOPES					OTOSCOPES & ACCESSORIES				VITAL SIGNS			5-Year Warranty**
	Wall Board	777 Wall Transformer	PanOptic Plus with Quick Eye Alignment	PanOptic Basic	Coaxial	MacroView Plus, Digitally-Capable	MacroView Basic	iExaminer® SmartBracket™ for MacroView Plus	Ear Spec Dispenser	Wall Aneroid	SureTemp Plus Thermometer	Connex® Spot Monitor/4400 Vital Signs Device*	
			118-3-US	118-2-US	11720-L	238-3	238-2	11840-IEX-US					
777 Integrated Wall Systems	777-PM3WCX-US	✓	✓	✓		✓		✓	✓			✓	✓
	777-PM3WAS-US	✓	✓	✓		✓		✓	✓	✓	✓		✓
	777-PM3WAX-US	✓	✓	✓		✓		✓	✓	✓			✓
	777-PM3WXS-US	✓	✓	✓		✓		✓	✓		✓		✓
	777-PM3WXX-US	✓	✓	✓		✓		✓	✓				✓
	777-PM2WCX-US	✓	✓		✓		✓		✓			✓	
	777-PM2WAX-US	✓	✓		✓		✓		✓	✓			
	777-PM2WXS-US	✓	✓		✓		✓		✓		✓		
	777-PM2WXX-US	✓	✓		✓		✓		✓				
	777-SM2WXX	✓	✓			✓	✓		✓				
Wall Transformer Sets	777-PM3XXX-US		✓			✓		✓					✓
	777-PM2XXX-US		✓				✓						
	777-SM2XXX		✓		✓		✓						

*Choosing these models will ensure the necessary mounting bracket ships with your wall board.
 Welch Allyn Connex Spot Monitor and/or Spot Vital Signs Device must be ordered separately.
 **5-year warranty is only applicable to the PanOptic Plus Ophthalmoscope and MacroView Plus Otoscope heads.

Welch Allyn® Green Series 777 Integrated Wall Systems and Transformers can help you stay efficient during your examinations. Select from wide-view and standard ophthalmoscopes and otoscopes, digital or manual blood pressure devices and integrated or stand-alone oral/axillary/rectal thermometry with Welch Allyn SureTemp® Plus Thermometer.



DIAGNOSTIC SETS

	OPHTHALMOSCOPES			OTOSCOPES & ACCESSORIES				POWER HANDLES (sets include one power handle)					5-Year Warranty**
	PanOptic Plus with Quick Eye Alignment	PanOptic Basic	Coaxial	MacroView Plus, Digitally-Capable	MacroView Basic	Traditional	iExaminer SmartBracket for MacroView Plus	Lithium-Ion Plus with built-in USB	Lithium-Ion	Nickel Cadmium	Hard Case	Soft Case	
	118-3-US	118-2-US	11720-L	238-3	238-2	250-2	11840-IEX-US						
71-PM3LXE-US	✓			✓			✓	✓			✓		✓
71-XM3LXE				✓			✓	✓			✓		✓
71-PX3LXERT-US	✓		✓					✓			✓		✓
71-PX3LXE-US	✓							✓			✓		✓
71-PM3LXES-US	✓			✓			✓	✓				✓	✓
71-PM2LXU-US		✓			✓				✓		✓		
71-SM2LXU			✓		✓				✓		✓		
71-SS2CXX			✓			✓				✓	✓		
71-SM2CXX			✓		✓					✓	✓		

**5-year warranty is only applicable to the PanOptic Plus Ophthalmoscope and MacroView Plus Otoscope heads.

Welch Allyn® Diagnostic Sets are the perfect solution to keep your exam tools with you for whenever you need them. Try our wide-view optics and fast charge power supplies built into our Plus models.



DESK SETS

All Welch Allyn Desk Sets shown below ship with a Universal Desk Charger and power supply.


Welch Allyn® Desk Sets are compact and perfect for smaller offices to support eye and ear examinations. Try our wide-view optics and fast charge power supplies built into our Plus models.



	OPHTHALMOSCOPES			OTOSCOPES & ACCESSORIES			POWER HANDLES (sets include two of the same handle)			
	PanOptic Plus with Quick Eye Alignment	PanOptic Basic	Coaxial	MacroView Plus, Digitally-Capable	MacroView Basic	iExaminer SmartBracket for MacroView Plus	Lithium-Ion Plus with built-in USB	Lithium-Ion	Nickel Cadmium	5-Year Warranty**
	118-3-US	118-2-US	11720-L	238-3	238-2	11840-IE-X-US				
71-PM3LDE-US	✓			✓		✓	✓			✓
71-PM2LDX-US		✓			✓			✓		
71-SM2LDX			✓		✓			✓		
71-SM2NDX			✓		✓				✓	
71-XXXLDE							✓			

**5-year warranty is only applicable to the PanOptic Plus Ophthalmoscope and MacroView Plus Otoscope heads.

Welch Allyn Ophthalmoscopes



	11720-L	118-2-US	118-3-US
Model	Traditional	PanOptic Basic	PanOptic Plus
Warranty	1 year	1 year	5 year
Larger field of view	—	20X	20X
Designed with DuraShock	—	Yes	Yes
Quick Eye alignment technology	—	—	Yes

Welch Allyn Oscopes



	250-2	238-2	238-3
Model	Traditional	MacroView Basic	MacroView Plus
Warranty	1 year	1 year	5 year
Larger field of view	—	3X	3X
Designed with DuraShock	—	Yes	Yes
Digitally-capable for image capture	—	—	Yes

Connect with your Hillrom representative today to discuss the best Welch Allyn Physical Exam & Diagnostic solution for your care setting.



hillrom.com