

NEEDLE COMPARISON CHARTS

RAZOR POINT™

PRECISION COSMETIC CONVENTIONAL CUTTING NEEDLES ▲

NEEDLE	CP MEDICAL	ETHICON	COVIDIEN	SHARPOINT	LOOK
13mm • 3/8 Circle	NPC-1	PC-1	DX-13	DSCM13	•
16mm • 3/8 Circle	NPC-3	PC-3	DX-16	DSCM16	•
19mm • 3/8 Circle	NPC-5	PC-5	DX-19	DSCM19	•



RAZOR POINT™


PRECISION REVERSE CUTTING NEEDLES ▼

NEEDLE	CP MEDICAL	ETHICON	COVIDIEN	SHARPOINT	LOOK
8mm • 3/8 Circle	NP-6	P-6	P-16	•	•
11mm • 3/8 Circle	NP-1	P-1	P-10	DSM11	C1
13mm • 3/8 Circle	NP-3	P-3	P-13	DSM13	C3
16mm • 3/8 Circle	NPS-3	PS-3	P-11	DSM16	•
19mm • 3/8 Circle	NPS-2	PS-2	P-12	DSM19	PC31
24mm • 3/8 Circle	NPS-1	PS-1	P-14	DSM24	PC33
26mm • 3/8 Circle	NPS	PS	•	•	•
8mm • 1/2 Circle	NP-2	P-2	•	•	•
16mm • 1/2 Circle	NPS-4	PS-4	P-24	•	PC30



REVERSE CUTTING NEEDLES ▼

NEEDLE	CP MEDICAL	ETHICON	COVIDIEN	SHARPOINT	LOOK
13mm • 3/8 Circle	NC-3	C-3	C-1	DS12	C17
16mm • 3/8 Circle	NFS-3	FS-3	C-12	DS16	C22
19mm • 3/8 Circle	NFS-2	FS-2	C-13	DS18	C6
24mm • 3/8 Circle	NFS-1	FS-1	C-14	DS24	C7
26mm • 3/8 Circle	NFS	FS	C-15	DS26	C8
30mm • 3/8 Circle	NFSL	FSL	C-16	DS30	C10
40mm • 3/8 Circle	NFSLX	FSLX	C-17	DS40	•
13mm • 1/2 Circle	NM-1	M-1	C-21	•	C26
15mm • 1/2 Circle	NX-8	X-8	•	•	•
17mm • 3/8 Circle	NJ-1	J-1	C-26	•	C20
23mm • 3/8 Circle	NX-1	X-1	C-23	HS24	•
26mm • 1/2 Circle	NCP-2	CP-2	GS-10	•	C21
37mm • 1/2 Circle	NCP-1	CP-1	GS-11	•	C12



[Needles are to scale]

NEEDLE COMPARISON CHARTS

STRAIGHT CUTTING NEEDLES ▽

NEEDLE	CP MEDICAL	ETHICON	COVIDIEN	SHARPOINT	LOOK
60mm	NKS	KS	SC-2	CS-1	•



TAPER POINT NEEDLES ○

NEEDLE	CP MEDICAL	ETHICON	COVIDIEN	SHARPOINT	LOOK
17mm • 3/8 Circle	NBB	BB	CV-15	•	•
17mm • 1/2 Circle	NRB-1	RB-1	CV-23	HR17	T28
22mm • 1/2 Circle	NCT-3	CT-3	GS-23	•	•
22mm • 1/2 Circle	NSH-1	SH-1	CV-25	HR22	T4
26mm • 1/2 Circle	NSH	SH	V-20	HR26	T
26mm • 1/2 Circle	NMO-6	MO-6	HGS-22	•	•
26mm • 1/2 Circle	NCT-2	CT-2	GS-22	•	T5
37mm • 1/2 Circle	NCT-1	CT-1	GS-21	•	T8
40mm • 1/2 Circle	NCT	CT	GS-24	•	T10



SOFTCUT NEEDLES ⊕

NEEDLE	CP MEDICAL	ETHICON	COVIDIEN	SHARPOINT	LOOK
17mm • 3/8 Circle	NV-4	V-4	KV-15	•	•
26mm • 1/2 Circle	NV-7	V-7	KV-7	•	•



CONVENTIONAL CUTTING STERNUM NEEDLES ▲

NEEDLE	CP MEDICAL	ETHICON	COVIDIEN	SHARPOINT	LOOK
48mm • 1/2 Circle	NCCS	CCS	SCC	•	•



[Needles are to scale]